

## Kontrollfacit PNV

ANALYS	Kontroll	-3SD	-2SD	MV	+2SD	+3SD	Lot
<b>Afinion/MicrosEmi</b>							
<b>CRPN</b> mg/L	PNV K4	9	12	16	21	23	K4 1909884, 2103819
	PNV K6	50	57	71	85	92	K6 1909885, 2103820
<b>HemoCue</b>							
<b>Hb</b> g/L	PNV HGL	73	75	78	81	83	HGL 13201, 14401
	PNV HGN	111	113	118	123	126	HGN 11902, 14402
<b>Accu-Chek inform II</b>							
<b>Glukos plasma</b> mmol/L	PNV CG1	1,8	2,0	2,4	2,8	3,0	CG1 10100331, 20103504
	PNV CG2	15,1	15,6	16,7	17,8	18,3	CG2 10100331, 20103504
<b>Clinitek Status</b>							
<b>Glukos</b>	PNV UA1			0			UA1 2101193, 2106072 UA2 2101194, 2106073
	PNV UA2		1		3		
<b>Ketoner</b>	PNV UA1			0			
	PNV UA2		1		4		
<b>Leukocyter</b>	PNV UA1			0			
	PNV UA2		1		4		
<b>Nitrit</b>	PNV UA1			0			
	PNV UA2			1			
<b>Protein</b>	PNV UA1			0			
	PNV UA2		1		3		
<b>Erythrocyter</b>	PNV UA1			0			
	PNV UA2		1		3		
<b>Förklaring:</b>	K4 Seronorm Liquid CRP-kontroll L-1 från Sero K6 Seronorm Liquid CRP-kontroll L-3 från Sero HGL, HGN Hemotrol från Eurotrol CG1, CG2 Accu-Chek Performa Control från Roche UA1, UA2 Liquid QC DipStrip Urinalysis Control från Triolab						